


# Adaptas' Next Generation MagneX<sup>®</sup> Micro Detector

Adaptas are proud to announce that we are bringing an all-new detector to market. Adaptas' next generation MagneX<sup>®</sup> Micro Detector encompasses award winning cross-field ETP MagneTOF<sup>®</sup> technology. This industry-leading technology provides the support needed for applications requiring a high performance miniature detector.

Some of the features of the new Micro detector include:

# High performance miniaturized particle/ion detection


Having the smaller size allows for improved accuracy and performance through miniaturized particle and ion detection. Our MagneX<sup>®</sup> Micro detector is approximately 1/3<sup>rd</sup> the size of detectors with comparable performance, and has higher linearity than detectors of a similar size.




# Cutting-edge reverse bias and magnetically permeable grid technology

The MagneX<sup>®</sup> Micro incorporates ETP's patented reverse bias technology to efficiently (>99%) focus electrons from the ion impact plate to the amplifying section via cycloidal electron transport [patent - WO2017015700A1].


#### Magnetically Permeable Grid






Transport of the electrons along the amplifying section is then optimised using a specialised ETP patented permeable grid [patent - WO2017132731A1].





# Superior dynamic range and gain for a compact detector

Capable of achieving 10  $\mu$ A linear output and 10<sup>6</sup> gain at 1400 V; the MagneX<sup>®</sup> Micro enables significantly increased dynamic range and exceptional gain characteristics for applications employing portable mass spectrometers such as residual gas and leak analysis, environmental and dedicated industrial sensors.



#### Fast pulse width

With a 1.5 ns FWHM single-ion pulse width; the MagneX<sup>®</sup> Micro can output ~4x the pulse height for the same gain (10 mV at 3x10<sup>6</sup>) of traditional pulse counting detectors. This extends its life and relaxes the gain requirement for pulse counting applications.

#### **Large Sensitive Area**

Detection efficiency has not been sacrificed. Despite its small form factor, the new detector features a

relatively large sensitive area. Utilizing advanced crossed-field optics, the MagneX<sup>®</sup> Micro transports and focusses electrons from a large ion-impact area to a miniaturized multiplication channel.

#### Ideal for established and emerging MS technologies

The MagneX<sup>®</sup> Micro is suitable for a wide range of MS applications. For the mini mass spec market, the Micro will decrease the size and increase the performance of next gen instrumentation and new innovations.



# **Specifications**

| Model number   | 14DM584                        |
|--|--------------------------------|
| Single-ion pulse width (FWHM)  | 1.5 ns                         |
| Multiple-ion pulse width (FWHM)  | ≤3.0 ns                        |
| Mechanical envelope size (nominal)   | 31 x 30 x 15.8 mm              |
| Input aperture diameter (nominal)  | 8.8 x 7 mm                     |
| Maximum sustained linear output current (typical)  | ~10 µA                         |
| Typical voltage for gain of 10 <sup>5</sup> (new multiplier)   | ~1200 V                        |
| Typical voltage for gain of 10 <sup>6</sup> (new multiplier)   | ~1400 V                        |
| Maximum multiplier –HV voltage (at end of life)<br>Note: Operating a new detector at maximum<br>voltage may result in damage to the detector | 3000 V                         |
| Maximum dark counts: -HV = -2000 V<br>(measured at 6E-6 Torr)  | <50 cts/min                    |
| Maximum dark current: -HV = -2000 V<br>(measured at 6E-6 Torr)   | <1 pA                          |
| Maximum operating pressure   | 10 <sup>-4</sup> Torr          |
| Long/short term storage requirements   | Protect from dust and humidity |
| External magnetic field strength   | ~200 G @ 5 mm (~20 G @ 15 mm)  |
| RoHS compliant   | Yes (exemption 7c-1)           |

### Dimensions

